GMM-H03A Body BMI Scale

Track your health in one scale

DESCRIPTION

- Body BMI scale utilizing ultrasonic to detect height, weight, and BIA*
- **BIA***: body & visceral fat, muscle percentage, bone mass, calories, and BMI
- Uses **software*** via PC to display body measurements
- Applicable for schools, hospitals, gyms, and clinics



H03A Body BMI Scale

det Volf-Mossaeren	alere a				< Edit Measurement M
Gom			×		s bait readin calent pr
	Bill Sod + Student Enter v D	D Search Skullers Date			
Granier + Next		D - D			
R Document r *	News	Determedited To			
& Dounload .r	Can List Sample CSF		Treads I		
Ritum /	T Student Exte Template Guide.pdf	23/7/2919 1044 am Ad	the Act		
Product in if	Citudent Data Temptata.cov	5/9/2021 11-21 are \$6	Setup		
and the second	e gere Cles List Sangle CIV			ram MOE Cockpit System	
	Regional Class Cat Sample CSV	Spen Cancel		ram MUE Cockpit System	
		agen Carca		Oute modified	
		· Marellar Usic	w	3/1/2021 12/21 em	
		Thatest Date	Weight 200 10330120946.COV Inreplate Guide pall	23/7.0019 10:44 am	
		01 Student Date	lemplate and	5/5/2021 11:01 are	
			Import		
dent: Hashh-Moraunament					- 0
dent Handle Mara norment	< End Session				- 0
-	< End Session	Stadent Data Inp			~ 0
Gmm	< End Session		et adreav v	Change Abset	- 0
-	< End Session		ednew v		- 0
-	< End Session	Name 1-4	sdrew ~ 5 Gerader N		- 0
-	< End Session	Name 14 Gaus P14 Height - G	adanır v 5 Gender N 6 Weight –	t - Kit	÷ a
Gmm	< End Session	Name 1-4 Gate P1-1	adanır v 5 Gender N 6 Weight –	t - Kit	- 8
-	< End Session	Name 14 Gaus P14 Height - G	adanır v 5 Gender N 6 Weight –	t - Kit	- 0
Gmm	< End Session	Name 14 Gaus P14 Height - G	adanır v 5 Gender N 6 Weight –	t - Kit	- 0
Gmm	< End Session	Name 14 Gaus P14 Height - G	adanır v 5 Gender N 6 Weight –	t - Kit	- 0
Gmm	< End Session	Name 14 Gaus P14 Height	adanır v 5 Gender N 6 Weight –	t - Kit	- 0
Gmm	< End Session	Name 14 Gaus P14 Height	adanır v 5 Gender N 6 Weight –	t - Kit	- 0
Gmm	< End Session	Name 14 Gaus P14 Height	adanır v 5 Gender N 6 Weight –	t - Kit	÷ 4
Gmm	< End Session	Name 14 Gaus P14 Height	adanır v 5 Gender N 6 Weight –	t - Kit	- 0
Gmm	< End Session	Name 14 Class Pj Height – co 10 523	odow v ; Gooder N 4 Weight – 65 Takes / N	1 - Kg 60	÷ 4
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- e
Gmm	< End Session	Andr	odow v ; Gooder N 4 Weight – 65 Takes / N	d on the	- 0
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	• 0
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- 4
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- 1
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- 4
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- 1
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- 4
Gmm	< End Session	Andr	staw v 5 Gender N 6 Weight - 65 Taken 2 N ew, please stan	d on the	- 1

Software via PC

FEATURES

- LCD display
- Weight measuring range: 5 200 kg (division: round off to nearest 50g)
- Height measuring range: 80 195 cm
- Power: Safety Mark AC/DC 3-pin power adaptor
- Customized software
- RS232 "scale-to-PC" communication cable
- Tempered glass weight base



SOFTWARE*

- Using RS232 cable as a communication medium between PC to scale, to run software
- Displays body measurements
- Data retrieved from software can be uploaded to microsoft excel for tracking



H03A base and LCD display



Visit testmeter.sg for more information!